

Tac-S101S BHK Serum-free Medium

Product Name: Tac-S101S

User Manual

Contents

Description	2
Application.....	2
Composition	2
Storage.....	2
Reconstitution of Powder Medium	2
Cryopreservation	3
Cell Recovery	3
Subculture Cells	4
Related Product	5

Description

Tac-S101S is a serum-free culture medium developed by Shanghai Biotechnology Co., Ltd. with independent intellectual property rights, targeting the growth and metabolism characteristics of BHK-21 cells. It is protein free and animal component free, suitable for high-density suspension culture of BHK cells, and supports efficient production of viruses such as pseudorabies, Newcastle disease, swine encephalitis B, and rabies.

Application

This product is intended for research or further manufacturing in the bio-manufacturing industry, but not for human or therapeutic use.

Composition

The IP rights of Tac-S101S medium formulation are owned by Shanghai BioEngine Sci-Tech Co., Ltd.

This medium contains:

- ☒ Carbohydrates, amino acids, vitamins, bulk salts, and trace elements.
- ☒ 8.7 g/L D-glucose, 1 g/L P188, 4 mM glutamine.

Not contain:

- ☒ Cytokines, antibiotics, HEPES and phenol red.
- ☒ Raw materials from animal sources.

Storage

- Store medium at 2-8°C, away from light.
- Once opened, the powder medium should be stored protected from moisture in a tightly sealed container.
- Do not use it after the expiration date or being damped.

Reconstitution of Powder Medium

Table 1 shows the preparation of Tac-S101S medium ^[1].

Ingredients	Concentration
Tac-S101S medium powder	24.72 g/L ^[2]
Sodium bicarbonate	2.44 g/L

Table 1. Preparation of Tac-S101S medium

- 1) Weigh 100% water of the final volume into the preparation container using pure water, ultrapure water, or water for injection at 20-30°C. Mix thoroughly without creating air bubbles.
- 2) Accurately weigh the corresponding mass of Tac-S101S medium at a concentration of 24.72 g/L and add it into the preparation container of 1) step. Stir well for 10-20 minutes.
- 3) Slowly adjust to pH 6.0-6.5 with 5-10 mol/L sodium hydroxide solution. The recommended amount of sodium hydroxide is 0.38 g/L. Stir for 10-20 minutes. At this point, the solution should be clear.
- 4) Weigh 2.44 g/L of sodium bicarbonate powder, add it slowly near the liquid level in the container, and stir for 10-20 minutes.

- 5) Adjust to pH 7.0-7.4 with sodium hydroxide or hydrochloric acid solution if the pH is beyond this range.
- 6) Pass the medium solution through a pore size of 0.22 or 0.2 μm sterile filter membrane, such as PES, using a pulse pump or compressed air (3-15 psi).
- 7) Use the prepared medium liquid immediately or store it in glass bottles, PET storage bottles, or disposable storage bags with an oxygen barrier membrane in a dark environment of 2-8°C. It's recommended for use within one month.

Note:

^[1] The above parameters (such as stirring time) are set for small-scale liquid preparation. Adjust these parameters for large-scale preparation based on container capacity to ensure full dissolution of dry powder.

^[2] The "g/L" unit denotes volumetric concentration (solute mass/solution volume).

Specifications of final liquid medium

Test	Unit	Specification
pH		7.0 – 7.4 ^[3]
Osmolality	mOsm/kg	310 – 350
Turbidity	NTU	< 1.00

Table 2. Specifications of final liquid medium

Note:

^[3] The pH buffer system of the product is carbon dioxide-sodium bicarbonate. The final pH value should be strictly controlled within the specific range outlined in Table 2. The following operations, such as prolonged reconstitution time or aeration in the bioreactor without pH control, can result in a gradual pH increase. There is a risk of metal ion precipitation when the pH value exceeds the upper limit.

Cryopreservation

- 1) Harvest cells in the mid-log phase of growth with >90% viability by centrifugation at 190×g for 5 minutes.
- 2) Prepare cryopreservation medium with 93% Tac-S101S medium and 7% DMSO on the day of use.
- 3) Resuspend cells in cryopreservation medium to a final viable cell density of 2.5-3.5×10⁷ cells/mL or as required.
- 4) Dispense aliquots of the cell suspension into cryovials.
- 5) Achieve cryopreservation in an automated or manual controlled rate freezing apparatus (0.5-1°C decrease per minute is suggested).
- 6) Transfer frozen cells to liquid nitrogen storage.

Cell Recovery

- 1) Rapidly thaw frozen cells in a 37°C water bath. Transfer to a super clean bench as soon as melted or with small ice crystals.
- 2) Transfer the vial content into a centrifuge tube containing 10 mL of prewarmed Tac-S101S medium. Harvest the cells by centrifugation at 190×g for 5 minutes and discard the supernatant.
- 3) Resuspend cells by prewarmed Tac-S101S medium to a viable cell density of 0.8-1.2×10⁶ cells/mL in a 125mL shake flask.
- 4) Incubate the shake flask at 37°C in a humidified atmosphere of 5% CO₂ in air on an orbital shaker platform rotating at 110-130 rpm (110 rpm for 50 mm amplitude; 130 rpm for 10 mm amplitude).
- 5) Cells should be sub cultured and adapted at least two passages. After the cell specific growth rate (or doubling time) reaches stability, subsequent operations can be carried out.

Subculture Cells

- 1) Ensure that the cell viability is >90%, and the growth rate is in mid-logarithmic phase prior to subculturing.
- 2) Calculate the volume of cell culture and prewarmed medium necessary to seed at $0.5-1.0 \times 10^6$ viable cells/mL in a shake flask.
- 3) Incubate at 37°C in a humidified atmosphere of 5% CO₂ in air on an orbital shaker platform rotating at 110-130 rpm (110 rpm for 50 mm amplitude; 130 rpm for 10 mm amplitude).
- 4) Subculture cells every two days according to the above steps.

Related Product

Product	Cat. No.	Form	Size	Packaging	Notes
Tac-S101S BHK Serum-free Medium	EXP0102003	Powder	10 L	Bag	- SF, PF, ADCF - Supports efficient vaccine production of pseudorabies, newcastle disease, Japanese encephalitis, rabies viruses, etc.
	EXP0102006	Powder	200 L	Barrel	
Tac-S201S BHK Serum-free Medium	EXP0119902	Powder	10 L	Bag	- SF, PF, ADCF - Support efficient vaccine production of rabies viurs
	EXP0119903	Powder	200 L	Bag	
LS206 BHK Serum-reduced Medium	EXP0112802	Powder	10 L	Bag	- SF, PF, ADCF - Lower serum addition(1%NBS) - Supports efficient replication of foot-and-mouth disease virus
	EXP0112803	Powder	200 L	Barrel	
SF201 BHK Serum-free Medium	EXP0101202	Powder	200 L	Barrel	- SF, PF, ADCF - Supports efficient production of the foot-and-mouth disease virus
	EXP0120002	Powder	10 L	Bag	
SF206 BHK Serum-free Feed Medium	EXP0105101	Powder	1 L	Bag	- SF, PF, ADCF - To be used with Tac-S101S or SF201 in fed-batch culture
	EXP0105104	Powder	20 L	Barrel	



Scan the QR code for more product information.

Stay tuned for more updates.

Tel: 86-21-68582660

Web: www.bio-engine.com.cn/EN

Shanghai BioEngine Sci-Tech Co., Ltd

Floor 4, Building 3, Lane 396, LvZhou Ring Road, Shanghai, China

Tel: 86-21-68582660

www.bio-engine.com.cn/EN



Shanghai BioEngine Sci-Tech Co., Ltd

Floor 5, Building 3, Lane 396, LvZhouHuan Road, Shanghai, China

Tel: 021-68582660

www.bio-engine.com.cn

